



A MEDICAL DEVICES MANUFACTURING COMPANY NEPHROLOGY | INTERVENTIONAL CARDIOLOGY | CRITICAL CARE

ISSUE 2.0 VOL- II : 2016-2017

## About the company

#### Introduction:

St.Stone Medical Devices Private Limited was incorporated in 2012 and start his planning for the one of the big project in North East India. Company has started his Official venture in 2014 by setup of Manufacturing Plant in North Eastern India. Althrought As a new born company having vast vision and mission to achieve as been granted the licence for manufacturing and Distribution of Medical Devices First time in North East India by CDSCO & State Drug Authority. We are into manufacturing & assembling of medical devices having vast range which falls under Interventional Cardiology Consumable, Nephrology, Aneasthesia and General Consumables.

#### Quality Accreditations:

We are Quality Certified ISO 13485:2012 and ISO 9001:2008 From UKAS Accreditation Body., CE Certification (Comp.), WHO-GMP Certification(Comp.) we are proud member of Pharmaexcil, All India Medical Device Association.

#### Research, Planning and Development:

As a growing company we don't have to work only on research and development. so we are working on research, planning and development., Our Hardcore team member are from Developed Companies from the same industry. who join us to make this happen that quality is the basic need in healthcare industry which help us to compete by making different strategies.

#### Manufacturing and Distributions:

As an growing company we need to think all the parameter which are essential for the Manufacturing of Medical Devices under WHO-TRS compliances. we although understand the major and minor necessities of Manufacturing and we step on it.

#### Manufacturing Area:

Our Manufacturing is as per standard of WHO Guidelines. Having 0.03 Um AHU Planted in Assembling, Manufacturing Area and Testing Lab. and 0.05Um AHU in Packaging, Sterilization and other areas of the premises.

#### Packaging Area:

Our Packing is as Classified into Class II as per WHO-TRS. Packaging Material is trans passes through Static Passboxes for much safety of the pouches and labels for not mixing and contaminations.

#### Sterilization:

St. Stone Having in house facility for Sterilization of Medical Device under EO Category. Sterilizer Calibration and Validation is been done in every intervals so we can assure on best quality for longer shelf life.

#### Distribution:

We are working on PAN INDIA so we can offer every one with our quality products. we are having International Distributor of ASEAN AND SAARC Counties

#### St. Stone Medical Vision:

St. Stone Medical Devices Vision is to provide Quality Healthcare at affordable to the people and making healthcare sector to be assured. St. Stone Medical Devices is not believing into making better Healthcare business, is believing into making Healthier Healthcare Industry for which St. Stone Medical Devices is working into.

#### St.Stone Medical Mission:

St. Stone Medical Devices is having aim to be the successful partner with many Healthcare Institution. The key to our success is our high degree of innovative drive coupled with the creativity of our highly qualified workforce

## EN ISO 13485:2012 | ISO 9001:2008 | ISO 14971:2007 | CE 1739 | WHO-GMP | DAC 1314QE28 | *Certified Company*

Member of AIMED | PHARMAEXCIL | FIEO I MANIIPUR CHAMBER OF COMMERCE & INDUSTRY

# DIALYSIS CONSUMABLE'S

ISSUE 2.0 VOL-II : 2016-2017

# **NEPHROLOGY / DIALYSIS CATHETER'S**

# **PLATICATH**<sup>TM</sup> HD Catheter

# **STRAIGHT/CURVE**

## Short Term Haemodialysis Catheters

An Optimum Compromise between competing Factors. Where a balance between:

- Flexibility & Firmness
- Softness and anti kinking

Double Lumen Catheter

- Material and durability
- Size and flow

#### Catheter Set Features

Long tapered tip dilators with stopper to avoid over insertion.

#### Set Contains:

- ▶ Polyurethane Catheter ▶ Vessel dilators ▶ Injection Cap
- ▶ Catheter Protector Pouch ▶ Patient Identification sticker

Anti-distorted and high flow rate can meet the requirement of femoral venous puncture.

# **PLATICATH™** DOUBLE LUMEN CATHETER

TYPE	TYPE DESCRIPTION		CURVE CODE
DOUBLE LUMEN	8.5 FR X 13 CM	PL8513S	-
DOUBLE LUMEN	8.5 FR X 13 CM		PL8513C
DOUBLE LUMEN	11.5 FR X 13 CM		PL11513C
DOUBLE LUMEN	11.5 FR X 15 CM	PL11515S	
DOUBLE LUMEN	11.5 FR X 16 CM	PL11516S	
DOUBLE LUMEN	11.5 FR X 20 CM	PL11520S	



#### Short Term Haemodialysis Catheters

- Flexibility & Firmness
- Softness and anti kinking
- Material and durability
- Size and flow

#### **Catheter Set Features**

Long tapered tip dilators with stopper to avoid over insertion.

#### Set Contains:

➢ Polyurethane Catheter ➢ Vessel dilators ➢ Injection Cap Catheter Protector Pouch 
Patient Identification sticker



# **PLATICATH™** TRIPLE LUMEN CATHETER

ТҮРЕ	DESCRIPTION	STRAIGHT CODE	CURVE CODE	
TRIPLE LUMEN	12 FR X 13 CM	-	PLT1213C	
TRIPLE LUMEN	12 FR X 15 CM	PLT1215S	-	
TRIPLE LUMEN	12 FR X 16 CM	PLT1216S	- (	6
TRIPLE LUMEN	12 FR X 20 CM	PLT1220S		

\*\*Above mention product are specified can be customized as per demand request\*

# **NEPHROLOGY / DIALYSIS CATHETER'S**

# **SHINECATH**<sup>TM</sup> Single Lumen HD Catheter



# SHINECATH<sup>™</sup> SINGLE LUMEN CATHETER

ТҮРЕ	DESCRIPTION	PRODUCT CODE
SINGLE LUMEN FEMORAL GREEN HUB	14GA / 13.5 CM	SHC14ANG
SINGLE LUMEN FEMORAL PINK HUB	14GA / 13.5 CM	SHC14AP

# Shinecath Single Lumen Femoral Catheter

An Optimum Compromise between competing Factors.

## Where a balance between:

- Flexibility & firmness
- Softness and anti kinking
- Material and durability
- Size and flow

## **Catheter Features:**

- ▶ Soft tip
- Small rotating fixation wing
- Extension lines stand repeated clamping and de-clamping
- Multiple threads female luer end
- Injection caps
- Variety of sizes, lengths and shapes

# **NEPHROLOGY / DIALYSIS CATHETER KIT**

# **PLATINUM<sup>+</sup>TM** Haemodialysis Procedure Tray



## Set Contains:

- ▶ Polyurethane HD Straight/Curve Catheter
- ▶ Guidewire with Advancer "J" Tip 70cm
- Catheter Protector Pouch
- ▶ 18Ga 7cm Introducer needle
- ▶ 5ml Guiding Syringe
- ▶ 5ml Disposable Syringe
- ▶ Vessel dilators ▶ Scalpel #11 ▶ Injection Cap
- ▶ Patient Identification sticker ▶ Drape Sheet
- ▶ Wound Dressing Pad






# **PLATINUM**+<sup>™</sup> DOUBLE/TRIPLE LUMEN CATHETER KIT

ТҮРЕ	DESCRIPTION	STRAIGHT CODE	CURVE CODE
DOUBLE LUMEN	8.5 FR X 13 CM	PP8513S	-
DOUBLE LUMEN	8.5 FR X 13 CM		PP8513C
DOUBLE LUMEN	11.5 FR X 13 CM		PP11513C
DOUBLE LUMEN	11.5 FR X 15 CM	PP11515S	
DOUBLE LUMEN	11.5 FR X 16 CM	PP11516S	
DOUBLE LUMEN	11.5 FR X 20 CM	PP11520S	
TRIPLE LUMEN	12 FR X 13 CM	-	PPT1213C
TRIPLE LUMEN	12 FR X 15 CM	PPT1215S	-
TRIPLE LUMEN	12 FR X 16 CM	PPT1216S	-
TRIPLE LUMEN	12 FR X 20 CM	PPT1220S	

# **NEPHROLOGY / DIALYSIS CONSUMABLE**



# SHINEWIRE SS Guidewire

## **SHINEWIRE™** ss guidewire

TYPE	DESCRIPTION	STRAIGHT CODE	3MM ANGLE CODE
GUIDEWIRE	0.021" X 50CM	SW2150ST	SW2150J
GUIDEWIRE	0.025" X 50CM	SW2550ST	SW2550J
GUIDEWIRE	0.032" X 50CM	SW3250ST	SW3250J
GUIDEWIRE	0.035" X 50CM	SW3550ST	SW3550J
GUIDEWIRE	0.038" X 50CM	SW3850ST	SW3850J
GUIDEWIRE	0.021" X 70CM	SW2170ST	SW2170J
GUIDEWIRE	0.025" X 70CM	SW2570ST	SW2570J
GUIDEWIRE	0.032" X 70CM	SW3270ST	SW3270J
GUIDEWIRE	0.035" X 70CM	SW3570ST	SW3570J
GUIDEWIRE	0.038" X 70CM	SW3870ST	SW3870J

Mercury Introducer Needle

Extremely low insertion forces

Smooth Guidewire Insertion

• Easy access and reduce tissue trauma.

Safety Cover for Unwanted Hazard
Variety of sizes, lengths and shapes.

## **Product Feature:**

- ▶ Accurate Measurement for Catheter Tracking ▶ Easy to insert and Removal
- ▶ Smooth Tip for insertion ▶ Angle and Straight Tip ▶ Kink Resistance
- Smooth Advancer for Dispensing of Guidewire Safety Cover for Unwanted Hazard Damage
- ▶ Variety of sizes, lengths and shapes.



# TYPEDESCRIPTIONPRODUCT CODECLEAR HUB-SQUARE18GA / 7CMMCCH1870PINK HUB-SQUARE18GA / 7CMMC1870PINK HUB-ROUND18GA / 7CMMCRH1870CLEAR HUB-SQUARE20GA / 3.81CMMC20381

# **NEPHROLOGY / DIALYSIS CONSUMABLE**





An arteriouvenus (AV) Fistula is an abnormal connection between an artery and a vein normally, blood flows from your arteries to your capillaries to your veins. Nutrients and oxygen in your blood travel from your cappillaries to tissues in your body. With an arteriovenous fistula, blood flows derictly from an artery into a vein by passing some capillaries.

Fistula needle are a crucial link between the patient and the dialysis machine. For this reason, all needles have to meet the requirements of highest quality, safety and comfort - both from the user's as well as from te patient's perspective.

# **STARFIST™** AV FISTULA NEEDLE

ТҮРЕ	DESCRIPTION	PRODUCT CODE
FIXED WING	16G/17G Lenght: Standard	SAVFXXXX
ROTATING WING	16G/17G Lenght: Standard	SAVRXXXX

\*\*\* XXX Stand for Abbreviation of Fr.Size and Lenght



An Optimum Compromise between competing Factors.

#### Where a balance between:

- ▶ Flexibility & firmness ▶ Softness and anti kinking.
- ▶ Material and durability ▶ Size and flow.

# **T-GUARD** TM Protector / Filter

#### Features:

- ▶ Efficient Filter ▶ PTFE Membrane ▶ 0.05Um Accuracy
- ▶ Smallest Size Packing ▶ Efficient Filter for Blood Transfusion


# **TGUARD™** TRANSDUCER PROTECTOR

ТҮРЕ	DESCRIPTION	PRODUCT CODE
0.5UM	PTFE	TG 03

\*\*Above mention product are specified can be customized as per demand request\*\*

# **CENTRAL VENOUS CATHETER KIT**

# **CENTVENT**<sup>TM</sup> CVC Catheter



The CVC is designed for placement in the Jugular Vein in critically ill patients. It's designed forlong term vascular access. For central veins catherization (subclavian, Jugular Veins) according to Seldinger method. Atraumatic metallic conductor with a J-shaped tip Prevents the vessel wall from injury.

# **CENTVENT™** CENTRAL VENOUS CATHETER KIT

ТҮРЕ	SIZES	LENGTH	PRODUCT CODE
	14G		CNTSXXYY
	16G	8CM, 13CM, 15CM,	CNTSXXYY
SINGLE LUMEN CATHETER	18G	16CM, 20CM	CNTSXXYY
	20G		CNTSXXYY
	22G		CNTSXXYY
	4.5 FR		CNTDXXYY
	5 FR	8CM, 13CM, 15CM, 16CM, 20CM	CNTDXXYY
DOUBLE LUMEN CATHETER	7 FR		CNTDXXYY
	8.5 FR		CNTDXXYY
	4.5 FR		CNTTXXYY
TRIPLE LUMEN CATHETER	5 FR	8CM, 13CM, 15CM, 16CM, 20CM	CNTTXXYY
	7 FR		CNTTXXYY

#### Set Contains:

- Central Venous Catheter
- SS Guidewire with Advancer "J" Tip 50cm
- Wound Dressing Pad
- ▶ 18Ga 7cm Introducer needle
- ▶ 5ml Guiding Syringe
- ▶ Vessel dilators ▶ Suture Wing
- ▶ Scalpel #11 ▶ Injection Cap
- Patient Identification sticker



\*\*\* XXYY Stand for Abbreviation of Fr.Size and Lenght

# **ANESTHESIA / SPINAL NEEDLE**

Spinal Introducer needles allow for easy insertion of the Spinal Needle and provide additional support.

Available in two types quincke and Pencil type Quincke Type needles provide a sharper bevel for low puncture force, while allowing for a tactile, reliable the precision-ground stylets are designed to minimized tissue coring and optimize gliding characteristics.



# MERCURY+<sup>™</sup> SPINAL NEEDLE

	ТҮРЕ	SIZES	LENGTH	PRODUCT CODE
		16G		SNXXYYQ
		18G		SNXXYYQ
		20G		SNXXYYQ
* *		22G		SNXXYYQ
T T	QUINCKE TYPE	23G		SNXXYYQ
		24G	40MM, 90MM, 120MM,	SNXXYYQ
		25G		SNXXYYQ
		26G		SNXXYYQ
		27G		SNXXYYQ
		16G	40MM, 90MM, 120MM,	SNXXYYQ
		18G		SNXXYYQ
		20G		SNXXYYQ
		22G		SNXXYYQ
	PENCIL TYPE	23G		SNXXYYQ
		24G		SNXXYYQ
		25G		SNXXYYQ
		26G		SNXXYYQ
		27G		SNXXYYQ

\*\*\* XXYY Stand for Abbreviation of NEEDLE GUAGE & LENGHT and Lenght

# CATHLAB CONSUMABLE'S

ISSUE 2.0 VOL-II: 2016-2017

www.stsmedicals.com



Manifolds are designed with easy-to-grip handles and clear polycarbonate bodies for visibility.

Medical's High & Medium Pressure Stopcocks and Manifolds are suitable for CathLab and other interventional applications and are designed to withstand 500 to 1200 psi pressure. Our products are available in 1 way and 3 way configurations with rotating High Pressure adapters or a standard swivel luer lock nut.

GANG	BODY	PRESSURE	ТҮРЕ	CATALOG CODE	
2 PORT	FULL BODY	600 PSI	<b>RIGHT ON</b>	STF225N	
2 PORT	FULL BODY	600 PSI	RIGHT OFF	STF225F	
3 PORT	FULL BODY	600 PSI	<b>RIGHT ON</b>	STF325NHB	
3 PORT	FULL BODY	600 PSI	RIGHT OFF	STF325FHB	
2 PORT	HALF BODY	500 PSI	<b>RIGHT ON</b>	STF20NBB	
3 PORT	HALF BODY	500 PSI	RIGHT ON	STF30NBB	

# STARFOLD<sup>™</sup> HIGH PRESSURE MANIFOLD 600PSI





## Y-Connector Push Click System

Y-Click Connector, provides an easy -touse hemostatic device for both diagnostic and interventional procedures. This ergonomically designed device significantly reduces the blood loss during interventional procedures. Among the Y-Click most unique



advantages are its self closing mechanism and semi -open position, as well as its one - handed operation.



## Y-Connector Push Pull System

Y-Push Pull Connector, which is used for coronary intervention, Push Pull head knob, which ensures minimal blood and precise navigation of devices in the vasculature. A specially engineered valve with micro vertical slits which facilitate smooth manipulation of device nullify blood loss.



## Y-Connector Touhy Borst System

Touhy Borst also know a Screw Type Y Connector, which is used for coronary intervention. It has a latest screw type head knob, which ensures minimal blood adjustable Thumb Wheel Mechanism pervents over tightening, loss of bloods and devices movement when turned. Allow for insertion of different size devices.

# **Y-STAR™ PTCA Y CONNECTOR**

ТҮРЕ	DESCRIPTION	PRODUCT CODE
Y-CLICK SET	Click Type Y Connector Hemostatis with extension line + Torque Device + insertion tool	YSTC9FI
PUSH PULL SET	Push Pull Y Connector Hemostatis with extension line + Torque Device + insertion tool	YSTPPS
TOUHY BORST	Screw Type Y Connector Hemostatis with extension line + Torque Device + insertion tool	STRTB9F





## **Inflation Device**

Uses while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon & deflate the balloon.

## **Product Feature:**

- ▶ High quality pressure sensing element makes Inflation Device pressure
- Pressure reading highly accurate even at extremes of pressure.
- Excellent graduation and error minimization for reading on grid and barrel
- Excellent eye to hand coordination by 30 ATM/Bar angled analouge pressure gauge
- Accurate pressure reading and high visibility.
- Feather Inflation Device ensures accurate repeatability of the measurements for several times.
- Lighter to handle and Easy to Operate

# **FEATHER™ INFLATION DEVICE**

	ТҮРЕ	DESCRIPTION	PRODUCT CODE
D	INFLATION DEVICE	Feather Inflation Device 30ml 30atm/bar with Stopcock	FID3030



## Pivot Introducer Sheath

#### Product Feature:

- ▶ Smooth transition of guidewire-to-dilator and dilator-to-sheath facilitates easy access
- Maintains uncompromised hemostasis
- Small outside diameter and large inside diameter reduces the size of the access site without compromising ease of device passage



# **PIVOT™** INTRODUCER SHEATH

ТҮРЕ	DESCRIPTION	PRODUCT CODE
FEMORAL SHEATH	5Fr to 8Fr / Length 7cm to 23cm	PIVFXXX
TRANSRADIAL SHEATH	5Fr to 8Fr / Length 7cm to 23cm	PIVTXXX

\*\*\* XXX Stand for Abbreviation of Fr.Size and Lenght



## Radial Compression Band

#### **Product Feature**

- Transparency ensures complete visual control from start to finish
- Dual compression balloons provide precise compression of the radial artery to ensure blood return-without compromising local nerve structure
- Innovative design makes it easy to use

# **STARBAND™** RADIAL COMPRESSION DEVICE


ТҮРЕ	DESCRIPTION	PRODUCT CODE
FEMORAL SHEATH	5Fr to 8Fr / Length 7cm to 23cm	PIVFXXX
FEMORAL SHEATH	5Fr to 8Fr / Length 7cm to 23cm	PIVTXXX



# **SHINEWIRE** TM PTFE Coated Guidewire

## Angiography PTFE Coated Guidewire Product Feature:

- Accurate Measurement for Catheter Tracking
- Easy to insert and Removal
- Smooth Tip for insertion
- ▶ Angle and Straight Tip
- ▶ Kink Resistance
- Smooth Advancer for Dispensing of Guidewire
- Safety Cover for Unwanted Hazard Damage
- ▶ Variety of sizes, lengths and shapes.



# **SHINEWIRE™** PTFE COATED GUIDEWIRE

OD SIZE	LENGTH	TIP	DESCRIPTION	CATALOG CODE
0.025″	150 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW25150J-SSMD
0.025″	150 CM	Straight "ST"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW25150ST-SSMD
0.032"	150 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW32150J-SSMD
0.032"	150 CM	Straight "ST"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW32150ST-SSMD
0.035″	150 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW35150J-SSMD
0.035″	150 CM	Straight "ST"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW35150ST-SSMD
0.038"	150 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW38150J-SSMD
0.038″	150 CM	Straight "ST"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW38150ST-SSMD
0.035″	260 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW32260J-SSMD
0.032"	260 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW35260J-SSMD
0.038″	260 CM	Angle "J"	3MM FIXED CORE, PTFE COATED GUDIEWIRE	SW38260J-SSMD

SHINEWIRE<sup>™</sup> is a range of complete guidewire, intended foe use of introduction spring winding, standard & special sizes, small diameter less grinding, specially tip forming, J-Tip, Straight Tip, Angle & other shape & depth marking.





## **Product Feature:**

- Small bore high pressure extension line.
- Suitable for high pressure monitoring
- Provided with male luer lock and female luer lock
- Luer lock compatible with all standard equipment.
- Luer lock connectors provide secure fitting.

# **LEICHIL<sup>TM</sup>** PRESSURE MONITORING LINE

ТҮРЕ	LENGTH	PRODUCT CODE
PRESSURE MONITORING LINE	50 CM	PML50-SSMD
PRESSURE MONITORING LINE	100 CM	PML100-SSMD
PRESSURE MONITORING LINE	150 CM	PML150-SSMD
PRESSURE MONITORING LINE	200 CM	PML200-SSMD

## Leichil<sup>™</sup> Pressure Monitoring Line 250-500psi

Pressure Monitoring Line is non expandable PVC Tubing with a ragid wall & non regid wall thickness which will ensure high quality wave forms with improved clarity. Available with both male-to-male & male-to-female connectors from 10cm to 200cms in length.




## Leichil<sup>™</sup> Extension Line with 3 way Stopcock

Extension Line with 3way stopcock is to be used with any of the Variot Infusion or transfusion devices in order to allow for attachment with any special or general purpose medical device. In order to improve ease of use also t provide an alternative channel for introduction of medicine to the patient.

<b>LEICHIL™</b> EXTENSION LINE WITH 3 WAY STOPCOCK			
ТҮРЕ	LEI	NGTH PROD	UCT CODE
EXTENSION LINE WITH 3 WAY S	ТОРСОСК 10	CM EXL10-S	SSMD
EXTENSION LINE WITH 3 WAY S	TOPCOCK 25	CM EXL25-S	SSMD
EXTENSION LINE WITH 3 WAY S	ТОРСОСК 50	CM EXL50-S	SSMD 🥢





# Leichil<sup>™</sup> High Pressure Monitoring Line 1000-1200psi

Pressure Monitoring Line is non expandable PVC Tubing with a ragid wall & non regid wall thickness which will ensure high quality wave forms with improved clarity. Available with both male-to-male & male-to-female connectors from 10cm to 200cms in length.



## Dose Control Syringe

## **Product Feature**

- > Transparent cylinder enhances visual effect of gauge mark and contrast level
- Special designed Silicone piston head and lubrication technology
- ▶ Easy to move plunger in cylinder and generate high pressure resistance
- Excellent eye to hand coordination
- Lighter to handle and Easy to Operate

# **EDHOU™** CORONARY CONTROL SYRINGE

ТҮРЕ	DESCRIPTION	PRODUCT CODE
MALE LUER LOCK	10 ML	EDU10
MALE LUER LOCK	12ML	EDU12

# COMPONENTS



# WE DO OEM FOR OUR SELECTED CUSTOMER AVAILABLE ALL CUSTOMIZED DISPOSABLE PRODUCT AS PER REQUEST

ISSUE 2.0 VOL-II : 2016-2017

Distribution & Represented By:-



Manufacturing Unit: **ST. STONE MEDICAL DEVICES PRIVATE LIMITED** G-03, Ground Floor, Dewlahland, Imphal, Manipur-India Licence No. : 01/DHSM/MFG (MD) Phone: 0385 2422827, 2422828 Email: info@stsmedicals.com

Head Office (Sales & Marketing): ST. STONE MEDICAL DEVICES PRIVATE LIMITED A-45/2, 2nd Floor, Devli Road, Khanpur, New Delhi-India Phone: 011 65002005, 8860554430 Email: stsmedicals@hotmail.com

ISSUE 2.0 VOL-II: 2016-2017